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GB English Instructions for Use GL-4M

Description


GL-4M is a programmable IR receiver with four relay outputs with alternating relay contacts. The receiver can be programmed with all GewaLink channels and one coded channel (4096). The receiver's relay can be programmed to function in a bistable or a monostable manner. The IR receiver has an internal IR detector that can be disconnected.

The following versions of the GL-4M are available:

4197	GL-4M	receiver for indoor use
419702	GL-4MU	receiver for outdoor use
419704	GL-4M OEM	receiver for mounting into equipment

Connections

All connections are attach by jacks, so that the connectors can be removed before connecting wires. The figure shows how wires are to be connected.

External Detector

The detector is the actual eye that receives the IR signal. A maximum of two external detectors can be connected at the same time. There is an internal detector on the circuit board of the receiver. This detector can be disconnected by cutting off the lines on the circuit board. This is shown in the figure.

The colours of the wires are to be connected as shown in the figure when connecting a detector. Make sure that the wires are correctly connected when connecting a detector.

Note: the range will be reduced if the detector eye is covered or hidden in some way, for example by curtains or by furniture.

Relay Connection

The receiver has four relays with alternating relay contacts. The maximum relay loading is specified in the Technical Specifications. Inductive loads, such as motors and relays, should be connected through a disturbance filter.

Power Supply

The receiver is powered with 12 - 24 VAC/DC (-10% - +20%).

Programming

The IR receiver can be programmed with GewaLink channels 0-127 (channel 62 has no function) or one coded channel (4096). There is one programming button for each relay. The selection is also made during programming whether a particular relay is to function in a monostable or a bistable manner.

Monostable functioning = The relay is activated as long as the key on the IR transmitter is held down.

Bistable functioning = The relay is activated when the key on the IR transmitter is held down and releases. The second time the key is held down and released, the relay will be deactivated.

Program the receiver as follows:

Press the selected button on the IR transmitter and press at the same time one of the programming buttons on the IR receiver. If the programming button is held down for approximately 1 second, the unit is programmed with a monostable function. If it is held down for approximately 2 seconds, the unit is programmed with a bistable function.


The receiver confirms the programming by activating the relay and causing the

indicator lamp to flash. Simply repeat the procedure if the selected channel or function is not satisfactory.

Warning! The relevant relay is activated during the programming. This may cause an inadvertent activation of some equipment.

Attaching the Lid

It is possible to fix the lid of the GL-4M model, with the included screws.

Cleaning

Dip a rag into water with a little washing-up liquid in it and wring it out well. This can then be used to clean the outer casing of the receiver.

Technical Specifications

Supply voltage:	12-24 VAC/DC (-10% - +20%)
Maximum current drawn when idle:	17 mA (3 detectors, relays inactive)
Maximum current drawn:	33 mA (3 detectors, all relays active)
Maximum relay load:	3A, 24 VAC/DC resistive load
Maximum number of detectors:	3, including internal detector
Programming channels:	GewaLink channels 0-127 (channel 62 is without function), one coded channel (4096) and Infracode.
Ambient temperature:	-20° - +35°C, (-4° - +95°F)
Surroundings:	Indoors: GL-4M, GL-4M OEM (IP-21) ” Outdoors: GL-4MU (IP-55)
Dimensions:	GL-4M: 96 x 96 x 36 mm ” GL-4M OEM 80 x 80 x 25 mm ” GL-4MU 100 x 115 x 57 mm

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Abilia AB, Kung Hans väg 3, Sollentuna, Box 92, 191 22 Sollentuna, Sweden