SE: Monteringsanvisning för IR-9MP Linak
GB: Installation manual for IR-9MP Linak

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GB: Installation manual for IR-9MP Linak

Important information

An electric power bed, wirelessly controlled by an IR transmitter, may cause personal injury if used or installed incorrectly.

The following precautions should be made:

Important information to the installation engineer

- Use only non-latching (monostable) control, do not use macro functions on an IR transmitter.
- If the user has physical difficulties at a certain angle, for instance difficulties in breathing, the minimum/maximum angle should be limited.
- On pneumatic control of an IR transmitter the sputum collector should be installed to collect superfluous sputum.
- Check that the user is able to reach the sip and puff mouthpiece in all angles on pneumatic control of the IR transmitter.
- The installing engineer shall give the following information to the user and assisting personnel.

Important information to users and assisting personnel

- On pneumatic control of an IR transmitter the sip and puff mouthpiece must be clean. Check that the mouthpiece will not be stopped up.
- IR transmitters are battery operated. Exchange battery at once when the battery indicates that it is running low or every 6 months.
- An IR receiver only works without problems if the IR signal from the IR transmitter has a clear view. Make sure nothing covers the eye of the IR receiver, for instance bedclothes.
- The IR receiver has no back up on power failure.


Description

IR-9MP LINAK is a programmable IR receiver with 8 relay outputs. A maximum 8 functions can be controlled on the bed, the head end, foot end, knee and the whole bed. The receivers are designed to be connected to the Linak control units CB-09, CB-10, CB-12 and Scanbed 750.

Connecting IR receiver

Connect the split cable to the control unit and the IR receiver and the manual control handle to the split cable. Hang the IR receiver on the wheel frame on the bed or on the headboard. Put the IR receiver in the direction where you get the best receiving conditions (not towards a wall). Inside the suspension device there is a Velcro to fasten the IR receiver in the right place.

Make sure no cables are squeezed! Fasten the redundant cabling with the enclosed strap. Check that the cables are long enough when the bed is in its highest position.


Programming channels

The default setting on the receiver is the channels 16-23 (see function table). The receiver is set to work with an IR transmitter with 8 functions. If another transmitter or other channels are to be used the IR receiver has to be re-programmed.

Do as follows to re-program:

Press a programming button on the circuit board on the IR receiver. Press, at the same time, the key on the transmitter that is to activate a function on the IR receiver. The IR receiver learns which channel is transmitted. If the programming button is pressed for approx. 1 sec. the function will be monostable. If the programming button is pressed for approx 4 sec. the function will be bistable. Note! When controlling a bed, the function must not be bistable. The indicator led on the receiver flashes to confirm the programming. Continue to program the next output etc.

Function table

Function	Programming button	IR-channel
Head up	1	16
Foot up	2	17
Knee up	3	18
Bed up	4	19
Head down	5	20
Foot down	6	21
Knee down	7	22
Bed down	8	23
X		X

Technical data

Supply voltage:	12-24 V DC
Current consumption standby:	11 mA
Max load/output:	140 mA from s.no: 497282 110 mA before s.no: 497282
Temperature range:	0° - +40°C
Programming channels:	GewaLink channels 0-61, 63-127 and 4096-code.

Compliance with the Medical Device Directive

The IR-9MP Linak is marketed as a medical device for persons with physical disabilities, and complies with Medical Devices Directive 93/42/EEC, which includes relevant EMC, electrical safety and risk analysis.

Abilia AB, Kung Hans väg 3, Sollentuna, Box 92, 191 22 Sollentuna, Sweden
Phone +46 (0)8-594 694 00 | Fax +46 (0)8-594 694 19 | info@abilia.se | www.abilia.se